

14 28.08  
**Signal** M-SERIES  
 Silicon Analyzer



**Silicon Analysis in Petroleum and Bio Fuels**

From gasoline to ethanol and toluene, the Signal bench-top analyzer delivers unprecedented precision and accuracy in quantitative analysis of silicon. The analyzer is based on XOS' MWD XRF technology platform (as applied in Sindie and Clara analyzers) ensuring a robust analysis solution for demanding petroleum and industrial environments.

**Application Areas:**

- Total silicon analysis in hydrocarbons and bio fuels.
- For use in refinery labs, pipeline terminals, additive plants, and inspection laboratories.

**Features and Benefits:**

- LOD: 0.5 ppm at 600 s.
- Dynamic Range: 0.5 ppm to 3000 ppm.
- Fits on any bench.
- Touch Screen user interface.
- User programmable measurement time: 30-900 s.
- No conversion gasses, heating elements, quartz tubes or columns.
- 75 W air-cooled excitation tube.

**Options:**

- LIMS compatible data output software.

**DASTEC S.R.L.**

**Distribuidores Exclusivos**

Buenos Aires, Argentina


Tel.: (54-11) 5352-2500

E mail: [info@dastecsrl.com.ar](mailto:info@dastecsrl.com.ar)


Web: [www.dastecsrl.com.ar](http://www.dastecsrl.com.ar)

# MWD XRF

**Monochromatic Wavelength Dispersive X-Ray Fluorescence (MWD XRF)** utilizes state-of-the-art focusing and monochromating optics to increase excitation intensity and dramatically improve signal-to-background over high power traditional WD XRF instruments. This enables significantly improved detection limits and precision and a reduced sensitivity to matrix effects. A monochromatic and focused primary beam excites the sample and secondary characteristic fluorescence x-rays are emitted from the sample. A second monochromating optic selects the silicon characteristic x-rays and directs these x-rays to the detector. MWD XRF is a direct measurement technique and does not require consumable gasses or sample conversion.



## Low Range Calibration



## Precision

Typical repeatability (r) and reproducibility (R) values in gasoline, at 95% confidence. 600 s measurement time.

Silicon Concentration (ppm)	r	R
2	0.4	0.7
5	0.5	0.8
8	0.6	1.0
15	0.8	1.4
100	2	4
500	5	10

## Product Specifications

<b>Dimensions</b>	37 cm (w) x 50 cm (d) x 34 cm (h)
<b>Power</b>	100-120 VAC, 47-63 HZ at 6.0 Amps/200-240 VAC, 47-63 HZ at 6.0 Amps
<b>Other Utilities</b>	Helium (10 psi maximum inlet pressure)
<b>Sample Cup Volume</b>	10 ml
<b>I/O Ports</b>	Ethernet 10/100 base T, RS232
<b>Optional Computer Interface</b>	Pentium, 100 MHz, 32 MB RAM/Windows 98 or newer operating system
<b>Ambient Temperature Requirements</b>	5-40°C (40-104°C)
<b>Dynamic Range</b>	Standard: 0.5 – 3000 ppm
<b>Measurement</b>	User selectable: 30-900 s
<b>Calibration</b>	8 calibration curves. Automatic and Manual Calibration functionality

**DASTEC S.R.L.**

**Distribuidores Exclusivos**

Buenos Aires, Argentina

Tel.: (54-11) 5352-2500

E mail: info@dastecsrl.com.ar

Web: www.dastecsrl.com.ar



**better analysis counts**

15 Tech Valley Drive • East Greenbush, New York 12061, USA • 518.880.1500 • Fax: 518.880.1510  
e-mail: info@xos.com • website: www.xos.com