



OPSIS LiquidLINE KjelROC Scrubber

Environment-friendly Digestion



KJELROC SCRUBBER

A Scrubber is recommended for fume control and in areas where water is expensive. The KjelROC Scrubber can save costs.

Reliable and Easy with a Small Footprint

The KjelROC Scrubbers are built to work in an environment with an excess of corrosive fumes and high heat. All materials are carefully selected.

- Cooler and containers in Borosilicate glass and connection tubes in Viton, in order to withstand the demanding environment and reduce maintenance costs
- Small footprint makes it easy to install the Scrubber close to the KjelROC Digestor
- Easy access to exchange containers, which are all made in high quality Borosilicate Glass.
- Accessible 6-step suction control to regulate suction power

NO WATER IS NEEDED

The KjelROC Scrubber does not require any tap water and there is a build-in pump for suction power.



Representantes / Distribuidores Exclusivos

Argentina Tel: (+54 11) 5352 2500 Email: info@dastecsrl.com.ar Web: www.dastecsrl.com.ar

🔚 Uruguay www.dastecsrl.com.uy

Paraguay www.dastecsrl.com.py

Opsis AB, Box 244 SE-244 02 Furulund Sweden Telephone +46 46 72 25 00 Telefax +46 46 72 25 01 E-mail info@opsis.se www.liquidline.se



All KjelROC Scrubbers can be supplied with optional configuration

- Water cooler, Water Accessory to even further increase the efficiency of the fume extraction. Can also be used with activated carbon filter.
- Water Trap for digestion of liquid samples





