

## Lab N<sub>2</sub> pressure meter **LAB.N2**



The newest ACM technology is the patented lab and at-line meter LAB.N2.

LAB.N2 determines non-destructively and precisely the container pressure of bottled still beverages with dosed N2 in PET and glass, like juices, waters and still wines. Longterm tests on the same sample are now possible. Also measurements on PET containers closed with sports caps are enabled.

Disposal of sample - common with invasive manometer methods - is economized.



New

## **Advantages**

No moving parts, maintenance-free Extremely fast, easy handling Loss-free sample treatment!

Sample loss-free and precise pressure determination now also for bottlers of non-sparkling beverages!



## Technical specifications

Measuring range **Accuracy** Repeatability **Measuring result Data interface** 

0-3,5 bar +/- 0,05 bar +/- 0.02 bar within I minute to ACM QUATROL



Representantes / Distribuidores Exclusivos

Argentina

Tel: (+54 II) 5352 2500 Email: info@dastecsrl.com.ar Web: www.dastecsrl.com.ar

Uruguay www.dastecsrl.com.uy

Paraguay www.dastecsrl.com.py



ACM GmbH Goldschlagstr. 172 1140 Wien / Austria

+43-1-7865866

**+43-1-7865866-20** 

office@acm.co.at

